VERI-Q

COVID-19 IgG/IgM Rapid Test

*MiCo BioMed*

## COVID-19 IgG /IgM

### COVID-19 IgG & IgM差分検出

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
*MiCo BioMed*

*Easy Answer to Life*

- The test detects COVID-19 antigens in direct nasal swabs or Nasopharyngeal Swab.

Test Window

C : Control

T : COVID-19 Antigen



Sample Well

CONFIDENTIAL

# PRINCIPLE OF THE ASSAY

MiCo BioMed

Easy Answer to Life

---

## Antigen Test Kit

---

Reaction Principle	Filter-through immuno-chromatography
Sample Type	Nasal Swab or Nasopharyngeal Swab
Test Results	Negative or Positive
Test Time	Less than 20 minutes for standard test
Storage condition	Cartridges : 2 ~ 30°C
	Extraction Reagent : 2 ~ 30°C

---

# PACKAGE CONTENTS

*MiCo BioMed*

*Easy Answer to Life*

- Individually pouched device with a desiccant (20 Test/kit)
- Extraction Reagent (5 mL X 2 bottles)
- Reagent Tube with Filter Cap (20 ea)
- Sterilized Nasal Swab (20 ea)
- Instructions for use (1 ea)

CONFIDENTIAL

# INTENDED USE

*MiCo BioMed*

*Easy Answer to Life*

The VERI-Q COVID-19 Ag Rapid Test is a lateral flow immunoassay intended for the qualitative **detection of nucleocapsid protein antigen from SARS-CoV-2 in direct nasal swabs** from individuals suspected of COVID-19 by their healthcare provider. Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen.

Antigen is generally detectable in nasal swabs during the acute phase of infection. Positive results **indicate the presence of viral antigens**, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses.

The VERI-Q COVID-19 Ag Rapid Test is **intended for use by medical professionals or trained operators** who are proficient in performing rapid lateral flow tests.


CONFIDENTIAL

# TEST PROCEDURE

*MiCo BioMed*

*Easy Answer to Life*

1. Add 300 uL extraction reagent (approximately 9 drops) into the reagent tube.
2. Add the sampled nasal/nasopharyngeal swab in the reagent tube with the extraction reagent.
3. Vigorously rotate and twist the swab against the side of the tube at least 10 times.
4. Squeeze the sides of the tube to obtain as much liquid as possible. Dispose of swab properly.
5. Apply the cap to the tubes with specimen, and hold the tube vertically, add 4 drops (about 100 uL) of the specimen without air bubbles into the sample well.
6. Start a timer. Read the result at 20 minutes. Test result should not be read after 30 minutes.

## Note:

For condensed samples, if the sample does not migrate to the membrane within 20 seconds, to apply one to two more drops of extract reagent to help the sample migrate on the membrane.

CONFIDENTIAL

# SPECIMEN PREPARATION


*MiCo BioMed*

*Easy Answer to Life*

Use the swabs supplied in the kit, other swabs may yield erroneous results. It is important to obtain as much secretion as possible. DO NOT USE FORCE while inserting the swab. The swab is designed to travel smoothly with minimum resistance.


## **Nasal Swab**

Carefully insert the sterile swab into the nostril that presents the most visible secretion or the nostril that is most congested. Using gentle rotation, push the swab until resistance is met at the level of the turbinates (less than one inch into the nostril). Rotate the swab a few times against the nasal wall and slowly remove from the nostril.

## **Nasopharyngeal Swab**

Carefully insert the sterile swab into the nostril that presents the most visible secretion or the nostril that is most congested. Keep the swab near the septum floor of the nose while gently pushing the swab with slow rotation into the posterior nasopharynx. Rotate the swab slowly several times and carefully withdraw.

## **Note:**

Specimens should be tested as soon as possible after collection. The sterile swabs should be opened immediately before use and specimens should be tested as soon as possible after collection.

All clinical specimens must be at room temperature before beginning the assay.

CONFIDENTIAL

# RESULT INTERPRETATIONS

MiCo BioMed

Easy Answer to Life

## Negative


A pink colored band appears only at the control line (C), indicating a negative result for COVID-19 Infections.

## Positive

A clear pink control band (C) and a detectable pink test band (T) appears, indicating a positive result for COVID-19 infections.

## Invalid

No visible band appears at the control line (C). Repeat with a new test kit. If the test still fails, please contact the distributor with the lot number.

The diagnostic accuracy of the VERI-Q COVID-19 Ag Rapid Test was evaluated for **clinical sensitivity and clinical specificity** by comparing results with RT-PCR, a standard diagnostic method.

## (1) Sample collection

- : 94 nasopharyngeal swab samples from SARS-CoV-2 infected patients
- : 201 non-SARS-CoV-2 infected nasopharyngeal swab samples

## (2) Test method

- : The test results were compared with RT-PCR analysis results.
- : SARS-CoV-2 Ag test were performed based on the Instruction for Use of VERI-Q COVID-19 Ag Rapid Test.

## (3) Test Results

- : VERI-Q COVID-19 Ag Rapid Test compared to RT-PCR results showed **97.3 % (95% CI: 86.2 % - 99.9 %) of diagnostic sensitivity and 100 % (95% CI: 97.3 % - 100 %) of diagnostic specificity in site 1.**  
And **85.7 % (95% CI: 73.8 % - 93.6 %) of diagnostic sensitivity and 100 % (95% CI: 94.7 % - 100 %) of diagnostic specificity in site 2.**

## A. Site 1 Evaluation Result.

Items		Reference method (RT-PCR)		
		Positive	Negative	Total
VERI-Q COVID-19 Ag Rapid Test	Positive	37	0	37
	Negative	1	133	134
	Total	38	133	171

- Diagnostic Sensitivity: 97.3 % (95% CI: 86.2 % - 99.9 %)
- Diagnostic Specificity: 100 % (95% CI: 97.26 % - 100 %)
- Positive predictive value (PPV): 100 %
- Negative predictive value (NPV): 99.3 % (95% CI: 95.1 % - 99.9 %)

## B. Site 1 Evaluation Result.

Items		Reference method (RT-PCR)		
		Positive	Negative	Total
VERI-Q COVID-19 Ag Rapid Test	Positive	48	0	48
	Negative	8	68	76
	Total	56	68	124

- Diagnostic Sensitivity: 85.7 % (95% CI: 73.8 % - 93.6 %)
- Diagnostic Specificity: 100 % (95% CI: 94.7 % - 100 %)
- Positive predictive value (PPV): 100 %
- Negative predictive value (NPV): 89.5 % (95% CI: 81.7 % - 94.2 %)

# PRECAUTIONS

*MiCo BioMed*

*Easy Answer to Life*

- For *in vitro* diagnostic use only.
- Do not re-use.
- Do not use if seal of the pouch is broken
- Do not use after the expiration date shown on the pouch.
- Do not mix and interchange different specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection while handling potentially infectious materials or performing the assay.
- Wash hands thoroughly after finishing the test.
- Do not eat, drink or smoke in the area where the specimens or kits are being handled.
- Clean up spills thoroughly with appropriate disinfectants.
- Handle all specimens as if they are potentially infectious.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and used devices in a proper bio-hazard container.
- Follow local, national or regional regulations when handling and disposal of the hazardous materials.
- Keep out of children's reach.
- Perform the test at 15 to 30 °C.
- Ensure that the test units are equilibrated at operation temperature before testing.

CONFIDENTIAL

# LIMITATIONS OF THE TEST

*MiCo BioMed*

*Easy Answer to Life*

- The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
- The contents of this kit are to be used for the qualitative detection of COVID-19 antigen from nasal swab, nasopharyngeal swab, nasal wash and nasal aspirate specimens.
- Neither the quantitative value nor the rate of SARS-CoV-2 antigen concentration can be determined by this qualitative test.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
- Failure to follow the Test Procedure and Interpretations of Test Results may adversely affect.
- Test Results must be evaluated in conjunction with other clinical data available to the physician.
- Negative test results do not rule-out possible other non-COVID-19 viral infections.
- Positive test results do not rule out co-infections with other pathogens.
- Positive and negative predictive values are highly dependent on prevalence. False negative test results are more likely during peak activity when prevalence of disease is high. False positive test results are more likely during periods of low influenza activity when prevalence is moderate to low.

CONFIDENTIAL