

Biocomma Limited

Pipet Tip Compatibility Chart(Continued)

Biocomma				TIP10—S	TIP10—LS	TIP20—S	TIP200—S	TIP200—LS	TIP1000—S	TIP1000—LS		
Pipettor Brand	Pipettor Make	Channels	Volume Range	0.1-10µL	0.1-10µL	1-20µL	1-200µL	1-200µL	100-1000µL	100-1000µL		
Eppendorf®	Reference	Single/multi-channel	0.1-2.5µL									
			0.1-10µL	●	●							
			0.5-10µL	●	●							
			2-20 µL “UltraMicro”									
			2-20µL			●	●	●				
			10-100µL			●	●	●				
			20-200µL					●	●			
			50-200µL					●	●			
	100-1000µL							●	●			
	Research	Single-channel	0.1-2.5µL									
			0.5-10µL	●	●							
			2-20µL			●	●	●				
			10-100µL			●	●	●				
			20-200µL				●	●				
			100-1000µL							●	●	
		Multi-channel	0.5-10µL	●	●							
			5-100µL			●	●	●				
			10-100µL			●	●	●				
			50-100µL				●	●				
			20-300µL				●	●				
			30-300µL					●	●			
	Research Pro	Single-channel	20-300µL					●	●			
			50-1000µL							●	●	
		Multi-channel	0.5-10µL	●	●							
			20-300µL				●	●				
	Research Plus	Single-channel	100-1000µL						●	●		
	Electronic	Single/multi-channel	0.5-10µL	●	●							
			5-100µL			●	●	●				
			2-200µL			●	●	●				
			50-1000µL							●	●	
	Single-channel	0.1-1µL	●	●								
		0.5-10µL	●	●								
		2-20 µL ‘Ultra Micro’										
		2-20µL			●	●	●					
		5-50µL			●	●	●					
		10-100µL			●	●	●					
		20-200µL				●	●					

BrandTech®	Transferpette®	Multi-channel	25-250µL				●	●			
			100-1000µL						●	●	
			0.5-10µL	●	●						
			2-20µL			●	●	●			
			2.5-25µL			●	●	●			
			5-50µL			●	●	●			
			10-100µL			●	●	●			
			20-200µL					●	●		
			30-300µL			●	●				
Thermo Scientific	Finnpipette™	Single-channel	0.2-2µL	●	●						
			0.3-3 µL	●	●						
			0.2-10µL	●	●						
			0.5-5µL	●	●						
			0.5-10µL	●	●						
			1-10µL	●	●						
			5-40µL			●	●	●			
			5-50µL			●	●	●			
			10-100µL			●	●	●			
			40-200µL				●	●			
			30-300µL				●	●			
			100-1000µL						●	●	
	200-1000µL						●	●			
		Multi-channel	1-10µL	●	●						
5-50µL				●	●	●					
50-300µL					●	●					
	Finnpipette Digital	Single/multi-channel	40-200µL				●	●			
	Finnpipette electronic	Multi-channel	1-10µL	●	●						
30-300µL					●	●					
100-1000µL							●	●			
			100-1200µL								
Gilson	Pipetman® Classic	Single-channel	P2, 0.1-2µL	●	●						
			P2, 0.2-2µL	●	●						
			P10, 0.2-10µL	●	●						
			P20, 2-20µL				●	●			
			P100, 10-100µL			●	●	●			
			P200, 20-200µL				●	●			
				P1000, 100-1000µL				●	●		
	Pipetman Neo	Single-channel	P2N, 0.2-2µL	●	●						
			P10N, 1-10µL	●	●						
			P20N, 2-20µL			●	●	●			
			P100N, 10-100µL			●	●	●			
P100N, 20-100µL						●	●				
P200N, 20-200µL				●	●						
Pipetman Neo	Multi-channel	P8x20N, 2-20 µL			●	●	●				
Pipetman Ultra	Single-channel	U200, 20-200 µL				●	●				
	Multi-channel	P200-M8, 20-200 µL				●	●				
Rainin	EDP3	Single-channel	E-MIC-10, 0.5-10 µL	●	●						
			E2-25, 2.5-25 µL			●	●	●			
			E2-100, 10-100 µL			●	●	●			
			E2-250, 25-250 µL				●	●			
		E2-1000, 100-1000 µL						●	●		
	Latch Mode	Single-channel	EP-M8-10, 0.5-10 µL	●	●						
			EP-M8-250, 10-250 µL			●	●	●			
			R2, 0.1-2 µL	●	●						
			R10, 0.2-10 µL	●	●		●	●			
			R20, 2-20 µL			●	●	●			
			R100, 10-100 µL			●	●	●			
			R200, 20-200 µL			●	●				

